2013 CATALOG Respiratory Protection

ACCESSORIES

Sperian is now Honeywell

Cleaning and Sanitizing

Refresher wipe pads provide a quick, easy method of removing perspiration and body oils from the respirator facepiece between scheduled cleanings.

These wipe pads make respirator usage more acceptable to the wearer by improving comfort. Shockwave™ is a disinfectant and cleaner that may be used to clean all reusable facepieces and their replacement parts. This concentrated liquid will make just the right amount of sanitizing solution for your needs when mixed with warm water. Each ounce will make one gallon of detergent that will sanitize, disinfect and deodorize surfaces.

SKU	Description	Packaging
7003	Refresher wipes with alcohol, 5" x 7"	Individually wrapped, 100 per box
7003L	Refresher wipes with alcohol, 8" x 10"	Individually wrapped, 100 per box
7003A	Refresher wipes, alcohol free, 5" x 7"	Individually wrapped, 100 per box
7003AL	Refresher wipes, alcohol free, 8" x 10"	Individually wrapped, 100 per box
80944A	Anti Fog wipes for North facepieces	Individually foil packaged, 10 per box
80995	Cleaner and disinfectant for respirator facepieces	10 oz. bottle
140082	Refresher wipes	Canister of 220 wipes
981803	Anti Fog wipes for Sperian facepieces	100 per box
981808	Anti Fog wipes, dry, for Sperian facepieces	100 per box



140082

Storage

OSHA and CSA require respirators to be stored in a clean and contaminant free environment. Protect your workers and your respirators by storing them properly.

SKU	Description	
79BAG	Nylon carrying bag North CFR-1, 7190 and mouthbit respirators	
77BAG	Nylon carrying bag half masks	
76BAG	Nylon carrying bag full facepieces	
200101	Tote Bag	
CA133	Molded carrying case for North Compact Air PAPR	
89986	Wall mount bracket for EBA	
83219	Wall Case, deluxe unit with shelf, for single SCBA	
940118	Zipper storage bag for half masks and full facepieces	











Call Customer Care at 800 430-5490.