

Detector Tube Pumps

Toximeter™ II Automatic Detector Tube Pump

The Toximeter™ II Automatic Detector Tube Pump makes the sampling process easier, allowing the user to preset the number of pump strokes from 1 to 250. Intrinsically safe, the automatic pump works with all MSA detector tubes. It can also be used as a sampling pump.

Toximeter™ II Automatic Detector Tube Pump

5142701T Toximeter™ II Automatic Detector Tube Pump



Kwik-Draw® Detector Tube Pumps

MSA's Kwik-Draw® and Kwik-Draw® Deluxe Pumps can be used with an assortment of MSA's detector tubes to spot-test atmospheres for a wide variety of toxic substances. They are designed for one-hand operation and consistent delivery of a sample draw volume of 100 milliliters (ml).

Kwik-Draw® Detector Tube Pumps

140904	Kwik-Draw® Basic Pump, with remote sampling adapter and carrying pouch
487500	Kwik-Draw® Deluxe Pump, with remote sampling adapter, carrying pouch and end-of-stroke indicator

Accessories

73067	Sampling line, 10 ft
73068	Sampling line, 25 ft
73069	Sampling line, 50 ft
87970	Remote Sampling Adapter, required for above sampling line
488780	Solvent-Resistant Sampling Line, 25 ft, with reel
488872	Tube Holder, required for 488780 line above
470321	Flue Gas Kit



Gas-Tester™ II H Detector Tube Pump

The Gas-Tester™ II H Pump is set for action by compressing the bellows. A pump stroke is started by pressing the release button. When the sample (100ml) is drawn through the tube, the end-of-stroke indicator changes color. An accurate measurement is obtained because the sample draw procedure itself is controlled only by the specifications of the pump and the flow resistance of the detector tube.

Gas-Tester™ II H Detector Tube Pump

5146760	Gas-Tester™ II H Detector Tube Pump
5146906	Gas -Tester™ II H Detector Tube Pump Service Kit



Detector Tube Handbook

Detector Tube Handbook

813929	Detector Tube Handbook
--------	------------------------